

INTERACT On MacOS

1 INTERACT On MacOS

The **Mangold INTERACT** Mac Edition has been tested on both **Intel**-chip, **M1**, **M2** and **M4** Chip devices, running **macOS 12, 13, 14*** and **15** (Monterey, Ventura, Sonoma, and Sequoia).

*) Requires an up-to-date version of **INTERACT**.

If you need to play multiple videos in sync, we recommend an M# Chip because it is far more powerful.

Note: *The macOS version is available for all INTERACT customers with a valid Scientific Service period.*

Licensing Method

We offer different licensing models, using either USB-dongles or personal license codes, and basically all work for both macOS and Windows.

If you have a **License Code** that is currently used on a Windows PC, contact us to move the license to a your mac.

Required Software Downloads

Download the required software through the following links:

➔ **License Software** - For both the **License Code** and the **License Key** method, the **CodeMeter Runtime kit** for macOS must be installed:
<https://www.wibu.com/support/user/user-software.html>

- Allow CodeMeter downloads.
- Mount the `CmRuntimeUser.dmg` by double-clicking the downloaded CodeMeter *.dmg file.
- Double-click the `CmInstall.pkg` file available on the disk image (dmg).

If you use your **CodeMeter License Key** on macOS for the FIRST time, check its compatibility at this point:

- Start the just installed **CodeMeter Control Center** app.
- Insert the license key.

If the **serial number of the stick pops-up** in the left pane of the Control Center, **your stick is compatible** so that you can continue with the installation of INTERACT.

If nothing appears, your stick is NOT compatible, at least not in its current state.

- Please get in touch with our [support team](#).

COMPATIBLE LICENSE KEYS

Customers that already own a **CodeMeter License Key** (dongle) need to check the serial number of the stick.

- Old sticks with a **serial number** starting with `1-.....` are NOT compatible with macOS.
You need to contact our [sales team](#) to discuss exchange conditions.
- Sticks starting with `2-.....` or `3-.....` might work right out of the box, if they are recognized by the CodeMeter Control Panel, which is part of the license communication software.

*If the **CodeMeter Runtime** kit on macOS does **NOT** recognize your license key, contact our [support team](#) to upgrade your license key during an online session on a Windows PC on which you have administrator rights.*

USING A LICENSE CODE

If you use a License Code, or your license key is visible in the Code Meter Control Center, you can install and run the **INTERACT** software on macOS:

➔ **INTERACT Software** - The **latest** macOS version of Mangold INTERACT is always available through these links, so make sure your service is up-to-date. Pick the version that fits your mac:


- INTERACT macOS **Intel**: <https://www.videoanalysis.net/setup/INTMac>
- INTERACT macOS **Silicon/ARM** (Mx):
<https://www.videoanalysis.net/setup/INTMacM>

IMPORTANT: Make sure your **license is still covered by an active service period** if you install the latest version!

- Download the INTERACT DMG file that fits your processor (Intel or ARM).

Installation on macOS

- Double-click the DMG file to mount it and drag the INTERACT app into the application folder:




- Start the INTERACT app, now available in your application folder.

The first time you start the application you'll most likely get this message:




- Click **Done**.
- Open the **System Settings** and navigate to the **Privacy & Security** section:
- Scroll down to the **Security** section where the message “INTERACT.app” was blocked.... is shown:



- Click **Open anyway**.
- If prompted, enter your credentials to confirm this action.
- After this, you may need to start INTERACT again, or it reopens automatically.

Anyway, you get a similar warning as the first time, except that you now have the option to click **Open Anyway**:

From now on, INTERACT will start directly – at least until the next update.




IMPORTANT: In addition to our INTERACT software, installing the license communication software [CodeMeter Runtime kit](#) for macOS from WIBU Systems is mandatory.

Python installation


Once you run INTERACT **for the first** time on your computer, **INTERACT** will most likely discover that Python 3.11 is missing and it will show the following prompt:



- Click OK.
- Double-click the `python-3.11.9-macos11.pkg`, if the installer is not started automatically.



Now the Python Installer wizard appears, prompting you to install this software.




- Proceed through the installer by clicking **Continue** a couple of times
- When you get prompted about the license agreement, click **Agree**.
- Wait until the installation is finished:



- Now restart the INTERACT application, so it can start downloading the required python packages. (This only works if Python is already available)
- Wait a few minutes so until all packages are downloaded.

You can verify this in the Terminal window that is open somewhere in the background:



```

pvmangold — 80x24
Installing pandas... Done!
Installing plotly... Done!
Installing pillow... Done!
Installing matplotlib... Done!
Installing seaborn... Done!
Installing holoviews... Done!
Installing pyvis... Done!
Installing wordcloud... Done!
Installing torch... Done!
Installing torchvision... Done!
Installing torchaudio... Done!
Installing openai-whisper... Done!
Installing GPUUtil... Done!
Installing pyannote.audio... Done!
Installing speechbrain... Done!
logout

Saving session...
...copying shared history...
...saving history...truncating history files...
...completed.

[Process completed]

```

- Wait until the [Process completed] message appears before you close this Terminal window!
- Now - if your license provides access to these features - all Python related routines, like Boxplots, Word Clouds and AI-based transcriptions should work as expected.

Troubleshooting

If Visualizations or the AI-based transcription do not work, although your license includes these features, there is probably a Python package conflict.

Especially if an older version of INTERACT (versions < 21.2) was previously installed on your Mac.

This is especially likely on Macs with a Silicon chip when earlier versions of INTERACT were installed, or whenever accidentally the wrong version was installed.*

* Also referred to as ARM. The ones with an M1, M2, M4 etc. processor

DELETE ALL OLD PYTHON COMPONENTS

- Delete the folder `Users/Shared/ Mangold Python`
- Delete the folder `Users/Shared/Mangold Python Data`
- Delete the folder `Users/Shared/ffmpeg`
- Delete the folder `HD/Library/Frameworks Python`
- Start **INTERACT** again and wait for all modules to download and install.



Mangold International GmbH
Graf-von-Deym Str. 5a
94424 Arnstorf, Germany

P +49 (0) 8723 / 978-330
F +49 (0) 8723 / 978-333

US / Canada: 404 537 2170

www.mangold-international.com